

QuikRead go ® Instrument software 8.2.2 release note

A new QuikRead go Instrument software version 8.2.2 will be released in August 2022.

The new software version has been developed and validated in compliance with IEC 62304:2006+AMD1:2015 Medical device software - Software life cycle processes and General Principles of Software Validation; Final Guidance for Industry and FDA Staff, Jan 11, 2002. Testing included measurement of 18 samples. All results passed.

An instrument with any of the previous software versions can be updated to this new software version. Version update does not delete old results from the memory of the instrument.

Features

Improvements

- GoDEX DT4x printer (Cat. no. 155129) and labels (Cat. no 155128) support added
- New TL-WN725N (Ver:3.0) WLAN adapter support added
 - o The Cat. no (148045) remains the same
- Fixed overlapping texts in iFOBT GUI layout
- Improved POCT1-A2 logging for customer support and R&D purposes
- SW 8.2.2 allows measurements on the date when the reagents expire

Availability of new software

- After release the new software version 8.2.2 is available on per customer basis from www.softwareupdate.quikread.com or can also be ordered as QuikRead go Software USB stick (Cat. no 137540).