

## QuikRead go Instrument software 8.3.0 release note

A new QuikRead go Instrument software version 8.3.0 will be released in February 2025.

The new software version has been developed and validated in compliance with IEC 62304:2006+AMD1:2015 Medical device software - Software life cycle processes and General Principles of Software Validation; Final Guidance for Industry and FDA Staff, Jan 11, 2002. Testing included measurement of 30 samples. All results passed.

An instrument with any of the previous software versions can be updated to this new software version. Version update does not delete old results from the memory of the instrument.

### Features

#### Improvements


- Fixed localized language translation errors in Chinese and Greek
- Fixed an issue where the Buf and Reag lots overlapped in result info screen
- Added support for special characters in SSID for wireless connection
- Added support for WPA-EAP networks
- Updated a hardware driver


#### Availability of new software


- After release the new software version 8.3.0 is available on per customer basis from [www.softwareupdate.quikread.com](http://www.softwareupdate.quikread.com) or can also be ordered as QuikRead go Software USB stick (Cat. no 137540).


U` à|â @âAâ [ { AD^AFI €FCBQ ] | Ç^âAX^|•â } AFAQ ] | Ç^âACB BâEG

**Signatures:**

<b>Document Review:</b>	I have reviewed the content of this document.	
Name:	<b>Kivirauma, Tapio</b>	
		2025-02-05 10:34:42 (UTC+00:00)
Electronically Signed in		Timestamp

<b>Document Review:</b>	I have reviewed the content of this document.	
Name:	<b>Dahal, Suman</b>	
		2025-02-05 10:39:06 (UTC+00:00)
Electronically Signed in		Timestamp

<b>Document Approval:</b>	I have Approved the content of this document and found no errors	
Name:	<b>Ramesh, Harish</b>	
		2025-02-05 12:02:06 (UTC+00:00)
Electronically Signed in		Timestamp

<b>Document Approval:</b>	I have Approved the content of this document and found no errors	
Name:	<b>Laukkonen, Jukka</b>	
		2025-02-05 14:54:01 (UTC+00:00)
Electronically Signed in		Timestamp

U` à|â @â|â| { AD|AF| €C|EQ| ] |! ç^âAX^!•â } AF|EQ| ] |! ç^â|CB| EB|EG